

Continue Learning with FirstHope

Complete Notes on the FirstHope App

Complete subject and semesters are available on FirstHope app only.

- Get complete subject and semester notes in the FirstHope app.
- Available for B.Pharm, D.Pharm, M.Pharm and Pharm.D students.

[Click here to get the FirstHope app](#)

Join the FirstHope Telegram Channel

Join FirstHope Notes and keep up with channel updates.

[Click here to join Telegram](#)

Instagram

@firsthope.official

[Click here to follow Instagram](#)

YouTube Channel

FirstHope-Pharmacy / Live streams

[Click here to watch the YouTube channel](#)

WhatsApp Channel

Follow FirstHope on WhatsApp.

[Click here to join the WhatsApp channel](#)

FirstHope Website

firsthopenotes.com

[Click here to visit the website](#)

Transverse Section and Powder Microscopy of Vinca Leaf

Experiment 4 - Practical Notes

Name: _____ Roll no.: _____ Date: _____

Aim

To examine a transverse section and powder of vinca leaf and record diagnostic features.

Principle

Authentication combines the arrangement of tissues in a section with diagnostic fragments in powder. *Catharanthus roseus* belongs to Apocynaceae. Material age, organ and processing can change the appearance; the complete pattern is more useful than one isolated cell.

Section Features

- Dorsiventral lamina with upper and lower epidermis, palisade tissue and spongy mesophyll.
- Midrib vascular region: distinguish xylem, phloem and surrounding ground tissue.
- Observe the stomatal pattern in an epidermal preparation. Anomocytic and anisocytic patterns have been reported; use the authenticated teaching specimen rather than imposing one universal pattern.

Powder Features

- Epidermal fragments with stomata.
- Palisade/mesophyll fragments and vessel elements from the veins.
- Other diagnostic inclusions actually demonstrated on the authentic slide.

Requirements

- Authenticated vinca leaf, labelled comparison slides and the prescribed microscope.
- Fresh razor used with a holder, slides, coverslips, water/glycerol mountant and instructor-prepared stains.

Procedure

1. Prepare a section

Use the instructor's safe cutting method to obtain a thin transverse section through the assigned organ. Transfer to water and select an intact section.

2. Mount

Mount in the assigned water/glycerol medium without folds or air bubbles. If staining is required, use the prepared reagent and prescribed rinse.

3. Low-power map

Observe the outline and tissue arrangement at low power. Draw the complete section with orientation and a clear title.

4. Diagnostic structures

Move to higher power to label the structures listed below. Record objective magnification; use stage calibration for any dimensions.

Pharmacognosy and Phytochemistry

5. Powder study

Mount a small amount of powder without a dense opaque layer. Scan several fields and draw at least three characteristic fragments.

6. Compare and conclude

Compare the section and powder with authenticated slides. Do not substitute a similarly shaped ornamental leaf for authenticated *C. roseus*. Record what was actually visible rather than copying every expected feature.

Observations and Record

Record	Entry
Botanical name / family / part	Catharanthus roseus / Apocynaceae / vinca leaf
Section, objective and stain	_____
Three powder fragments and magnification	_____
Diagnostic agreement / limitations	_____

Worked Example

The numbers below are illustrative; enter your own observations in the record.

- Magnification alone does not give a fragment length. If calibration is 4 μm per eyepiece division and a fragment spans 12 divisions, its length is 48 μm .

Result

State whether the observed section and powder support the stated botanical material, naming the features seen and any unresolved difference.

Precautions

- Use authenticated, labelled teaching material and instructor-prepared reagents. Wear the protection specified for the chemical method.
- Compare with a known control and record actual observations. A single colour, precipitate or odour description cannot authenticate a crude drug.
- Do not taste botanical material or directly inhale vapours; dispose of solvents and staining reagents in their assigned waste containers.
- Use teacher-prepared staining/clearing reagents. Do not make chloral-hydrate or corrosive staining mixtures independently.

Viva Questions and Answers

1. Why study both section and powder?

The section preserves tissue arrangement, while powder reveals fragments retained after grinding.

2. Why examine several fields?

A single field may miss a diagnostic feature or overrepresent a contaminant.

3. How should an uncertain feature be recorded?

Describe it as uncertain and request comparison; do not force an identification.

Continue Learning with FirstHope

Complete Notes on the FirstHope App

Complete subject and semesters are available on FirstHope app only.

- Get complete subject and semester notes in the FirstHope app.
- Available for B.Pharm, D.Pharm, M.Pharm and Pharm.D students.

[Click here to get the FirstHope app](#)

Join the FirstHope Telegram Channel

Join FirstHope Notes and keep up with channel updates.

[Click here to join Telegram](#)

Instagram

@firsthope.official

[Click here to follow Instagram](#)

YouTube Channel

FirstHope-Pharmacy / Live streams

[Click here to watch the YouTube channel](#)

WhatsApp Channel

Follow FirstHope on WhatsApp.

[Click here to join the WhatsApp channel](#)

FirstHope Website

firsthopenotes.com

[Click here to visit the website](#)