

Complexometric titration

Complexometric Titration

Introduction:

- Complexometric titration is a volumetric analysis technique used to determine the concentration of metal ions by forming a complex with a chelating agent (titrant).
- It is commonly applied in fields such as chemistry, environmental science, and pharmaceutical analysis, particularly for quantifying metal ions, especially transition metals.

Principle:

- The titration is based on the formation of stable complexes between metal ions and chelating agents.
- Chelating agents, such as EDTA, have multiple donor atoms (e.g., oxygen, nitrogen) that form coordinate bonds with metal ions, creating complex ions.

Key Components:

1. Analyte:

- The metal ion (or ions) whose concentration is being determined.

2. Titrant (Complexing Agent):

- A chelating agent added in a known volume and concentration to form complexes with the metal ions.
- Common agents include EDTA and DTPA.

3. Indicator:

- Signals the endpoint, often through a color change.
- Examples include Eriochrome Black T or Murexide, which change color when the metal ion forms a complex with the chelating agent.

Applications:

- Determination of metal ion concentrations in environmental samples (water, soil).
- Analysis of trace metals in pharmaceuticals.
- Quality control in food and beverage production.
- Monitoring heavy metal pollution in industrial effluents.
- Research in analytical chemistry and metallurgy.

Methods for Endpoint Detection:

1. Visual Indicators:

- Indicators like Eriochrome Black T change color when they form a complex with metal ions.
- A second color change indicates the endpoint when all metal ions are complexed.

2. Potentiometric Detection:

- Measures potential difference between electrodes.
- The endpoint is marked by a sharp change in potential as all metal ions react with the chelating agent.

3. Conductometric Detection:

- Monitors the electrical conductivity of the solution.
- The endpoint is reached when conductivity changes significantly as metal complexes form.

4. pH-Metric Detection:

- Tracks pH changes during titration.
- A sudden pH shift indicates the endpoint when all metal ions are complexed.

PM Indicator in Complexometric Titration

- The term **PM indicator** refers to the **Patton-Reeder indicator**, used in complexometric titrations.
- It is not a category, but a specific indicator named after its developers, **Patton and Reeder**.

Common Indicators Used in Complexometric Titrations:

1. Eriochrome Black T (EBT)
2. Murexide (Ammonium Purpurate)
3. Calcon (1-hydroxy-3-(m-sulfonatophenylazo)-2-naphthoic acid)
4. Solochrome Dark Blue (SDB)
5. Xylenol Orange

Classification of Complexometric Titrations

1. Direct Titration:

- In this method, metal ions in the sample are titrated directly with a chelating agent, such as EDTA.
- The endpoint is detected using an indicator like Eriochrome Black T or by instrumental methods (e.g., potentiometric or conductometric detection).
- **Direct titrations** are the simplest and most commonly used complexometric titrations.

2. Back Titration:

- Back titration is used when direct titration is challenging, such as when the endpoint is not well-defined, or the reaction is slow.
- In this method, the metal ions are first reacted with an excess of the chelating agent.
- The unreacted chelating agent is then titrated with a standard solution of a metal ion (e.g., zinc or magnesium) of known concentration.
- The endpoint is detected using the same methods as in direct titration.

3. Replacement (Displacement) Titration:

- Used when the metal ions do not readily form a stable complex with the chelating agent or when the endpoint is difficult to detect.
- A metal ion (such as calcium or magnesium) that readily forms a stable complex with the chelating agent is added to displace the analyte metal ions.
- The displaced metal ions are then titrated with the chelating agent, and the endpoint is detected similarly to direct titration.

4. Alkalimetric Titration:

- Although not a type of complexometric titration, **alkalimetric titration** refers to titrations used to determine the acidity or alkalinity of a solution using a strong acid or base.
- Also known as **acid-base titrations**, these are not typically categorized under complexometric titrations.

This arrangement clarifies the different types of complexometric titrations, highlighting the appropriate classification and method of detection for each.

Metal ion indicators

- **Metal ion indicators**, also known as **chelating indicators** or **metallochromic indicators**, are organic compounds that form colored complexes with metal ions.
- They are used in **complexometric titrations** to detect the endpoint, which occurs when all the metal ions in the sample have reacted with the chelating agent.
- Indicators are selected based on their ability to form **weak complexes** with the metal ions being analyzed.

How They Work:

- The **formation of a colored complex** between the metal ion and the indicator causes a visible color change.
- During the titration, the **chelating agent** (e.g., EDTA) has a stronger affinity for the metal ions and displaces the indicator, breaking the weak complex and causing a color change.

- When all metal ions have reacted with the chelating agent, any **excess chelating agent** reacts with the indicator, causing another color change that signals the **endpoint** of the titration.

This concise explanation covers the role and function of metal ion indicators in complexometric

Masking and demasking reagents

Introduction

- Masking and demasking are techniques used to selectively protect or reveal certain ions or functional groups during chemical reactions or analysis.
- These techniques are particularly useful in complex mixtures to control the reactivity of specific components and prevent unwanted side reactions.

A. Masking

- Masking involves temporarily protecting a specific ion or group by forming a stable, less reactive complex with a masking agent.
- This allows selective reactions to occur without interference from the masked species.
- **Example:** In complexometric titrations, a masking agent can bind certain metal ions, preventing them from reacting with the titrant while allowing other ions to be titrated.

B. Demasking

- Demasking is the reverse process where the masked ion or group is released by altering conditions (e.g., pH) or adding a demasking agent.
- **Example:** After completing the titration of unmasked ions, demasking agents can be added to release the previously masked ions, allowing them to be titrated separately.

Masking and Demasking Reagents

Common Masking Reagents

- **Ethylenediaminetetraacetic acid (EDTA):** Masks many metal ions by forming stable complexes.
- **Cyanide ions (CN⁻):** Masks metals like silver (Ag⁺) and copper (Cu²⁺) by forming cyanide complexes.
- **Dimethylglyoxime (DMG):** Masks nickel (Ni²⁺) by forming a stable red complex.
- **1,10-Phenanthroline:** Masks iron (Fe²⁺) by forming a stable complex.

Common Demasking Reagents

- **Hydrogen peroxide (H₂O₂):** Demasks chromium (Cr³⁺) from its EDTA complex by oxidation.
- **Oxalic acid:** Demasks calcium (Ca²⁺) from fluoride by forming calcium oxalate.
- **Ammonium hydroxide (NH₄OH):** Demasks aluminium (Al³⁺) and iron (Fe³⁺) from fluoride complexes by raising the pH.



FREE PREVIEW

[Click here to visit our website](#)

PHARMACY PREMIUM NOTES

Complete subject
and semesters are
available on
FirstHope app only.



[Scan or click here](#)

Click here to get the app



Available on Google Play

JOIN THE FIRSTHOPE COMMUNITY

Join our Telegram channel

Click here to join Telegram



[Scan or click here](#)

@FirstHopeNotes

Instagram

@firsthope.official

[Click here to follow >](#)

YouTube channel

FirstHope-Pharmacy / Live streams

[Click here to watch >](#)

WhatsApp channel

Follow FirstHope on WhatsApp

[Click here to join >](#)

FirstHope website

firsthopenotes.com

[Click here to visit >](#)

B.Pharm / D.Pharm / M.Pharm / Pharm.D

Preview notes. Complete subjects and semesters in the FirstHope app.