

Ka Ke Half-life Cmax AUC AUMC and MRT

Experiment 7 - Practical Notes

Name: _____ Roll no.: _____ Date: _____

Aim

To analyse the assigned ka ke half-life cmax auc aumc and mrt experiment or teaching dataset.

Principle

A concentration-time profile provides observed and model-derived parameters. Estimating absorption/elimination constants requires an appropriate model and sufficient sampling.

Requirements

- Institution-assigned experimental method and coded in-vitro or supplied pharmacokinetic data.
- Calibrated analytical readings or supplied concentration table, calculator and plotting tools.

Procedure

1. Plot data

Plot concentration versus time and terminal log concentration versus time; identify observed Cmax and its sampling time.

2. Estimate Ke and half-life

Fit an appropriate terminal region and calculate Ke from its natural-log slope.

3. Assess Ka

For the assigned one-compartment first-order oral model, use the prescribed residual or fitting method; do not call the terminal slope Ka without justification.

4. Calculate areas

Use the assigned integration method for AUC and AUMC, showing measured and extrapolated components.

5. Report MRT

Divide AUMC by AUC on the same integration basis and state the route/model interpretation.

Observations and Record

Record component	Value
Observed Cmax / Tmax	_____
Terminal slope / Ke	_____
Ka and model	_____
AUC / AUMC	_____
Half-life / MRT	_____

Biopharmaceutics and Pharmacokinetics

Worked Example

The numbers below are illustrative; enter your own observations in the record.

- $K_e = -\text{slope of } \ln(C) \text{ versus } t$; $t_{1/2} = \ln(2)/K_e$. For terminal $\log_{10}(C)$, $K_e = -2.303 \times \text{slope}$.
- Linear-trapezoidal AUC segment = $(C_i + C_{i+1})/2 \times (t_{i+1} - t_i)$. AUMC applies the same stated integration approach to t_C .
- For monoexponential terminal extrapolation: $AUC_{\text{tail}} = C_{\text{last}}/K_e$; $AUMC_{\text{tail}} = t_{\text{last}} C_{\text{last}}/K_e + C_{\text{last}}/K_e^2$. $MRT = AUMC/AUC$.
- For the applicable first-order one-compartment oral model: $T_{\text{max}} = \ln(K_a/K_e)/(K_a - K_e)$, when K_a and K_e differ.
- Examples: $K_e = 0.2 \text{ h}^{-1}$ gives $t_{1/2} = 3.466 \text{ h}$. $AUC = 40 \text{ mg}\cdot\text{h/L}$ and $AUMC = 120 \text{ mg}\cdot\text{h}^2/\text{L}$ give $MRT = 3 \text{ h}$. A segment from 0 h/0 mg/L to 1 h/4 mg/L gives $2 \text{ mg}\cdot\text{h/L}$.

Result

Dataset / method: _____. Calculated outcome and units: _____. Comparison: _____. Assumption or limitation: _____.

Precautions

- Record units, sampling times, dilution and the assumptions of the selected model.
- Use supplied human/animal data or the approved in-vitro teaching method; no new human/animal dosing or surgical preparation is prescribed.
- Keep an illustrative calculation separate from the actual result and do not use it as an individual patient dosing instruction.

Viva Questions and Answers

1. Why state model assumptions?

The derived parameter is valid only when the data and model satisfy them.

2. What accompanies a pharmacokinetic result?

Its units, dataset, calculation method and relevant limitation.

Continue Learning with FirstHope

Complete Notes on the FirstHope App

Complete subject and semesters are available on FirstHope app only.

- Get complete subject and semester notes in the FirstHope app.
- Available for B.Pharm, D.Pharm, M.Pharm and Pharm.D students.

[Click here to get the FirstHope app](#)

Join the FirstHope Telegram Channel

Join FirstHope Notes and keep up with channel updates.

[Click here to join Telegram](#)

Instagram

@firsthope.official

[Click here to follow Instagram](#)

YouTube Channel

FirstHope-Pharmacy / Live streams

[Click here to watch the YouTube channel](#)

WhatsApp Channel

Follow FirstHope on WhatsApp.

[Click here to join the WhatsApp channel](#)

FirstHope Website

firsthopenotes.com

[Click here to visit the website](#)